

allard_{INT}

-Peek-a-boo!



ToeOFF® The original



Now – pediatric eggs in our basket!

kiddieGAIT™

The newest addition to the ToeOFF Family of Products
Floor reaction control and dynamic assist for toddlers



The newest addition to the ToeOFF® Family of Products!

Pediatric gait poses special challenges due to the dynamic muscle tone that often presents in much of this population. Not only are there the usual ground reaction forces involved in gait biomechanics, but also top-down spasticity driven forces exacerbate gait deviations and deformities.

Now, KiddieGAIT offers innovative options that have never been available for AFO management of these challenges. Functional environments can be created that supplement gait function instead of immobilizing and inhibiting that function.

Start off on the right foot – give toddlers an orthosis that will allow their little feet to move in a more fluid and natural biomechanical gait pattern.

- Carbon composite with anterior shell, lateral strut, and an engineered footplate with open calcaneus.
- Designed to allow integration of your custom foot orthotic to help control ground-up forces.
- Lateral strut can also serve as a base to affix a T-Strap to control excessive eversion at the ankle.
- Anterior shell to assist in management of proximal deficits by helping to overcome either knee hyper-extension or crouch gait.

The ToeOFF® carbon composite AFO has proven itself invaluable for more than ten years in helping manage adult gait deficits, including those with some tone involvement. Now you can extend many of the benefits of this design concept to manage the lower extremity functional challenges of this younger population.

INDICATIONS

- Footdrop
- Gait deviation secondary to proprioceptive deficit (Either unstable or low-tone gait)
- Toe-walker with no midfoot collapse
- Low Tone Crouch Gait
- Spina Bifida
- Cerebral Palsy
- Muscular Dystrophy
- MMC

CONTRAINDICATIONS

- Lacking ROM towards dorsiflexion (need at least 5° dorsiflexion past neutral)
- Very rigid foot structure
- Quadriceps spasticity
- Fixed postural Genu Valgum or Genu Varum
- Fixed postural Pes Valgus or Pes Varus

LIMITATIONS

- Knee hyperextension



Give them the stability they need without locking up their little ankles!

Benefits

- Functional heel-to-toe gait when used in conjunction with appropriate foot orthotic.
- Anterior support provides opportunity to enhance proprioceptive response.
- Lightweight—so important for these little feet.
- Can use less brace to do the same work as traditional plastic AFOs.
- Footplate gives spring assist - KiddieGAIT gives back energy.
- Anterior support adds stability and opportunity to manage knee extension/flexion moment.
- Softer and more dynamic heel rocker, which creates softer loading to reduce knee flexion moment.
- Therapy friendly - able to use SMO only for therapy, or with KiddieGAIT when stability and/or dynamic assist is required.
- Wearer compliance is increased - lightweight functional support allows more fluid and productive gait (less tiring).

FOOT INTERVENTIONS

KiddieGAIT should always be combined with an additional orthotic, designed to control the position of the foot. It is important that the foot position is as close to neutral as possible when bearing weight. The most common foot related problems like pronation, supination, pes varus, pes valgus should be corrected with this additional orthotic.

When spasticity is present, it is generally recommended that this orthotic include a deep heel cup to further encourage the heel-to-toe gait.

The following are guidelines only. When the foot related problems are:

Mild = use a Foot Orthotic that is deep and firm

Moderate = UCB (Deep Foot Orthotic with trim line below ankle)

Severe and/or tight achilles tendon = SMO (Supra Malleolus Orthosis or Molded Inner Boot)

For low tone (hypotonia) and pronation, consider a soft SMO, such as SureStep™.

SHOE

CONSIDERATIONS

- When spasticity is present, a rocker sole with approximately a 10mm heel height. For hypertonia, consider a flat sole but build up the heel - either under the footplate or externally to accommodate the heel rise built into KiddieGAIT.
- Removable insole
- Firm counter
- Secure closure



Floor Reaction Control and Dynamic Assist

Combine the best of a solid ankle AFO with the best of an articulating AFO. KiddieGAIT offers movement dynamics similar to an articulating, while offering resistance similar to a solid. The carbon graphite shell adds stability, while its lightweight and energy reflecting properties encourage more dynamic gait.





Functional Heel-to-Toe Gait!

kiddieGAIT™

KiddieGAIT offers a carbon composite framework for you to build an orthosis to encourage functional heel-to-toe gait. While an SMO (molded inner boot) or UCB controls the foot, the KiddieGAIT floor reaction frame adds dynamic response, lower leg anterior support and stability, and/or opportunity for management of knee extension/flexion moment. This combination offers a unique opportunity for you to customize function to meet specific patient needs.





Allard USA expresses deepest gratitude and appreciation to those who helped to make this brochure possible: All the beautiful children, their parents and siblings who participated in "We Kicked Up Our Heels" event and allowed us to share their photos; Orthotic Specialists Inc.; Boulder Mountain Therapy, Inc.; Photographer John Michael Tarchala



kiddieGAIT™

When spasticity is present, KiddieGAIT should always be fitted with a semi-rigid to rigid foot orthotic, UCB, SMO (Molded Inner Boot), that provides proper positioning and control of the foot. For low tone (hypotonia) and pronation, consider a soft SMO, such as SureStep™.

SureStep. Supra Malleolus Orthosis. Foot Orthotic

Item no	Product	Size	Height	Footplate Length
28372	KiddieGAIT	Small*	220mm 8 5/8"	160mm 6 1/4"
28372	KiddieGAIT	Medium*	257mm 10 1/8"	180mm 7 1/16"
28372	KiddieGAIT	Large*	295mm 11 5/8"	200mm 7 7/8"
28397	SoftKIT	Small	Patent no.: ToeOFF; ToeOFF Short; ToeOFF Fantasy; ToeOFF NFR; BlueRocker; BlueRocker NFR; KiddieGAIT; KiddieGAIT NFR; U36950; BE1005297; BE1114626; DK1005297; DK1114626; FI1005297; FI1114626; FR1005297; FR1114626; IE1005297; IE1114626; IT1005297; IT1114626; CA2279225; CNZL97181689.1; NL1005297; NL1114626; NO313656; PL194247; CH1005297; CH1114626; ES1005297; ES1114626; GB1005297; GB1114626; SE1005297; SE1114626; DE69709184.8-08 DE69732541.5-08; US5897515; ATE210417; ATE289187 Ypsilon; Ypsilon NFR GE60208889.5-08; GB1379201; IT1379201; US6887213	
28397	SoftKIT	Medium		
28397	SoftKIT	Large		

*Specify left or right