

A New development within the ToeOFF® Family!

NAVIGAIT™ is a new alternative for users with a mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the user's specific needs.

- Developed to better assist user's with a mild foot drop to gain better dorsiflexion
- The padded shell is positioned comfortably on the shin.
- Removable interface.
- The amount of dorsiflexion is adjusted through the position in the lacing and with the velcro attachment on the shin plate.
- Low profile design makes it invisible under slacks.
- NAVIGAIT™ is available in two sizes.
- The plastic shin part can be trimmed with scissors if needed.

NAVIGAIT™

Intended use and indications

NAVIGAIT™ is intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple sclerosis, Post polio syndrome.

Contraindications

NAVIGAIT™ should not be used by patients who present with spasticity or severe oedema.

Limitations

NAVIGAIT™ does not address malpositions of the foot or foot deformities. These should be cared for with shoes and/or orthotics.

For other indication, we refer you to the other products within the ToeOFF®-Family.






One important part in our product development is to listen to the market and collect ideas and requirements that are given to us from users and professional within the orthopedic industry.

The development of NAVIGAIT™ is a unique cooperation between users, CPO and Allard International. As one step closer to our slogan "Support for better life", the continued development of the product range is important as we strive to fulfill the needs of individuals with physical challenges.

Peter Allard
CEO

Item No.	Description	Size	Height
28410 0011	NAVIGAIT™	Small/Medium	27,5 cm
28410 0013	NAVIGAIT™	Large/XLarge	31,5 cm
28412 0011	Extra strap w/coil	Small/Medium	
28412 0013	Extra strap w/coil	Large/XLarge	

patent pending



September 2010 ©