

Selection® Radialis Finger and thumb parts



Individual support for hand and fingers – in combination with Selection® wrist orthoses

Selection® Radialis Finger and Thumb Parts are designed to provide targeted support for the hand and fingers when used together with Selection® wrist orthoses: 35204, 35205, 35209, 35222, 35223, and 35305.

They can be used for patients with Radial Nerve Palsy, where nerve damage affects grip function. In such cases, they are recommended in combination with the longer MultiOpen wrist orthosis, which offers additional wrist support during rehabilitation.



Available in smaller sizes for children

The Thumb Part is also available in XS and XXS, making it compatible with our children's wrist orthosis. This makes it particularly well-suited for children with conditions such as cerebral palsy, who often require more tailored solutions.

Radialis Finger Part

The Finger Part consists of two sections, making it easy to put on and adjust for an optimal fit. It attaches securely to the soft Velcro surface on the dorsal side of the wrist orthosis and provides comfortable support to the fingers.

Radialis Thumb part

The Thumb Part is individually adjustable – unlike orthoses with an integrated thumb section and offers full grip support throughout the rehabilitation phase.

Flexible and individualized solutions

The Finger and Thumb Parts are available separately and can be used individually or combined, depending on the patient's needs.

Key benefits of Selection® Radialis:


- Low profile and high comfort – consistent with the Selection® wrist orthosis range
- Flexible adjustment for individual needs
- Suitable for both children and adults
- Ideal for Radial Nerve Palsy and other impairments affecting hand and finger function



Selection® Radialis Finger

Art. No.	Size	L/R	Measurement
35217 1011	Small	Left	17 - 20 cm
35217 1012	Medium	Left	20 - 23 cm
35217 1013	Large	Left	23 - 27 cm
35217 2011	Small	Right	17 - 20 cm
35217 2012	Medium	Right	20 - 23 cm
35217 2013	Large	Right	23 - 27 cm


Measurement:
Circumference MCP 2-5



Selection® Radialis Thumb

Art.no.	Size	L/R	Measurement
352161009	XX-Small	Left	1,5-2 cm
352161010	X-Small	Left	2-3 cm
352161011	Small	Left	3 - 4 cm
352161013	Large	Left	5 - 6 cm
352162009	XX-Small	Right	1,5-2 cm
352162010	X-Small	Right	2-3 cm
352162011	Small	Right	3 - 4 cm
352162013	Large	Right	5 - 6 cm

Measurement: Length from joint center MCP1 to joint center IP1.





Support for Better Life

Everyone should be able to live their life to the fullest, regardless of their mobility challenges. With innovative solutions developed in close collaboration with healthcare professionals and patients, we strive to provide Support for Better Life.

allard

Allard International

Karbingatan 38, 254 67 Helsingborg, Sweden | +46 42 25 27 00 | info@allardint.com | allardint.com