

VISSION textile orthosis

Now we are launching a new unique material, developed for a new series thumb and wrist assortments. The soft comfortable inside of cotton has been combined with a durable surface of polyester, nylon and lycra.

Design and pattern development has focused on the anatomical fit to provide smooth and consistent orthoses.

Vission Wrist Orthosis Short

A slightly shorter wrist orthosis, made of sturdy and durable resilient material with a soft and comfortable inside. The model is completely open and easy to wear. The orthosis has a good anatomical fit inside the hand, over the back of the hand and around the muscles of the thumb. This model is without any volar spoon. It has an open adjustable thumb strap. The model is 13 cm long on the volar side.

Intended use: Supports and relieves the wrist after injury, surgery, abrasion, agitation, rheumatic problems, chronic pain conditions, joint wear or over strain. Can also be used for prevention in heavy or monotonous work.



| Item no Left | Item no Right | Colour | Sizes | Wrist circumference |
|--------------|---------------|--------|----------|---------------------|
| 46030 1009 | 46030 2009 | Black | XX-Small | 11 - 12 cm |
| 46030 1010 | 46030 2010 | Black | X-Small | 12 - 13 cm |
| 46030 1011 | 46030 2011 | Black | Small | 13 - 14 cm |
| 46030 1012 | 46030 2012 | Black | Medium | 14 - 16 cm |
| 46030 1013 | 46030 2013 | Black | Large | 16 - 19 cm |
| 46030 1014 | 46030 2014 | Black | X-Large | 19 - 22 cm |

Vission Thumb Orthosis with Mouldable Plastic



This thumb orthosis has a plastic insert made from low temperature thermoplastic. It should be shaped to fit the patient's anatomy by dipping into a water bath (72°C). It is made of sturdy and durable resilient material with a soft and comfortable inside. The model is completely open and easy to wear. The orthosis has a good anatomical fit inside the hand, over the back of the hand and around the muscles of the thumb. It provides a stable support. It has an open and adjustable thumb strap. The model is 18 cm long on the thumb side.

Intended use:

This product is intended to restrict movements in CMC and MCP I joints and to position the thumb in a functional position by an individually moulded plastic thumb support. Mild to Moderate CMC and MCP I instability, Ligament Injuries, Morbus de Quervain, Tendinitis, Degenerative diseases, Trauma.

| Item no Left | Item no Right | Colour | Size | Wrist circumf. |
|--------------|---------------|--------|---------|----------------|
| 46024 1010 | 46024 2010 | Black | X-Small | 12 - 13 cm |
| 46024 1011 | 46024 2011 | Black | Small | 13 - 14 cm |
| 46024 1012 | 46024 2012 | Black | Medium | 14 - 16 cm |
| 46024 1013 | 46024 2013 | Black | Large | 16 - 19 cm |
| 46024 1014 | 46024 2014 | Black | X-Large | 19 - 22 cm |

Low Temperature
Thermoplastic Insert



Vission Thumb Orthosis with stay

This thumb orthosis have a plastic stay on the thumb side to support the thumb and the surrounding structures. The material under the stay gives good padding to avoid pressure on the joint or muscles. The material is durable and stable, with a comfortable inside. The thumb strap is open for adjustability. The short model is 13 cm long and the longer model is 16 cm on the thumb side.

Intended use:

This product is intended to be used for instability in the thumb joints (CMC- or MCP I joints) and to position the thumb in a functional positions. Indications such as inflammation, joint stress and for rheumatic degenerations fo the short model and also Trigger Finger, Morbus De Quervain for the long model.



SHORT

| Item no Left | Item No Right | Model | Colour | Size | Wrist circumference |
|--------------|---------------|-------|--------|----------|---------------------|
| 46028 1009 | 46028 2009 | Short | Black | XX-Small | 11 - 12 cm |
| 46028 1010 | 46028 2010 | Short | Black | X-Small | 12 - 13 cm |
| 46028 1011 | 46028 2011 | Short | Black | Small | 13 - 14 cm |
| 46028 1012 | 46028 2012 | Short | Black | Medium | 14 - 16 cm |
| 46028 1013 | 46028 2013 | Short | Black | Large | 16 - 19 cm |
| 46028 1014 | 46028 2014 | Short | Black | X-Large | 19 - 22 cm |
| 46023 1009 | 46023 2009 | Long | Black | XX-Small | 11 - 12 cm |
| 46023 1010 | 46023 2010 | Long | Black | X-Small | 12 - 13 cm |
| 46023 1011 | 46023 2011 | Long | Black | Small | 13 - 14 cm |
| 46023 1012 | 46023 2012 | Long | Black | Medium | 14 - 16 cm |
| 46023 1013 | 46023 2013 | Long | Black | Large | 16 - 19 cm |
| 46023 1014 | 46023 2014 | Long | Black | X-Large | 19 - 22 cm |



LONG

Vission Wrist Orthosis with Stay



A wrist orthosis in standard model made of sturdy and durable resilient material with a soft and comfortable inside. The model is completely open and easy to wear. The orthosis has a good anatomical fit inside the hand, over the back of the hand and around the muscles of the thumb. It provides a stable support and has a moldable plastic spoon on the volar side. It has an open adjustable thumb strap. The model is 19 cm long on the volar side.

Intended use:

Intended for when the wrist needs support and off-loading in cases such as inflammations, work related strain or when rheumatic and/or degenerative changes are present.

| Item no Left | Item no Right | Colour | Size | Wrist circumference |
|--------------|---------------|--------|----------|---------------------|
| 46029 1009 | 46029 2009 | Black | XX-Small | 11 - 12 cm |
| 46029 1010 | 46029 2010 | Black | X-Small | 12 - 13 cm |
| 46029 1011 | 46029 2011 | Black | Small | 13 - 14 cm |
| 46029 1012 | 46029 2012 | Black | Medium | 14 - 16 cm |
| 46029 1013 | 46029 2013 | Black | Large | 16 - 19 cm |
| 46029 1014 | 46029 2014 | Black | X-Large | 19 - 22 cm |

Vission Wrist Orthosis with Thumb

This model have the same basic pattern as the standard wrist, and also a support around the thumb. A soft stay on the thumb side gives a light support. Both the wrist spoon and thumb stay can be removed. It has an open adjustable thumb strap. The model is 19 cm long on the volar side.

Intended use:

Intended to restrict movements in CMC and MCP1 joints and to position the thumb in a functional position and to reduce flexion in the wrist. Indications such as Mild to Moderate CMC and MCP1 instability, Ligament Injuries, Tendinitis, Degenerative diseases, Trauma.

| Item no Left | Item no Right | Colour | Size | Wrist circumference |
|--------------|---------------|--------|----------|---------------------|
| 46031 1009 | 46031 2009 | Black | XX-Small | 11 - 12 cm |
| 46031 1010 | 46031 2010 | Black | X-Small | 12 - 13 cm |
| 46031 1011 | 46031 2011 | Black | Small | 13 - 14 cm |
| 46031 1012 | 46031 2012 | Black | Medium | 14 - 16 cm |
| 46031 1013 | 46031 2013 | Black | Large | 16 - 19 cm |
| 46031 1014 | 46031 2014 | Black | X-Large | 19 - 22 cm |

