


Product information

# ***ELEMENTS***

*for Retailers*





# GENERAL INFORMATION

ELEMENTS is a line of products designed and developed to meet the daily life of the ordinary person that suffers from pain in joints and muscles. These braces are only for one person use.

The materials in ELEMENTS have been carefully chosen and tested to be functional during long term use - during the day without giving to much warmth.

The foam in the ELEMENT-material is made from fully cured polyurethane. The textile has been heat laminated towards the foam. No harmful adhesives have been used during the production process.

The packaging is made from polypropylene. This material can (depending on local recommendations) normally be handled in the household waste.

The polyester material in ELEMENTS is Öko-Tex Approved and certified for class II, standard 100. This means that it is safe to use directly on skin.

Always instruct the customers not to use body lotions or oils on the skin when using the brace.

If your customer has no previous experience using a brace we recommend to instruct the patient to start slowly with a few hours per day and then increase the use gradually.

Please encourage your customer to read the instructions supplied with the product carefully.

Good Luck!  
Allard Int.

## *ELEMENTS*


# ANKLE SUPPORT

## Instructions for ELEMENTS ANKLE Item No: 34100

**Intended use:** This product is intended for treatment of mild and temporary conditions such as ankle pain and ankle sprains. The support can be used with or without the elastic straps. It gives a light compression and keeps the ankle joint warm. The movement of the ankle joint is limited by the elastic straps.

**Sizes:**

34100-0011	Small	20 - 23 cm
34100-0012	Medium	23 - 26 cm
34100-0013	Large	26 - 29 cm



**Material:**

Beige/Black outer material:	Polyamid and Spandex
Foam:	Polyurethane
Velcro:	Nylon
Elastic straps:	Polyester, Nylon and Lycra
Edging:	Nylon and Spandex
Decoration:	Nylon and Spandex

**Washing Instruction:**

Close all Velcro prior to laundering. Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

**Contraindications:**

Severe edema, severe swelling, neuropathy, foot problems due to diabetes, leg ulcers.



# ELEMENTS


# ELBOW SUPPORT

## Instructions for ELEMENTS ELBOW SUPPORT Item No: 34101

**Intended use:** This product is intended for treatment of mild and temporary conditions such as pain and inflammation in the elbow. The support gives a light compression and keeps the elbow warm.

**Sizes:**

34101-0011	Small	23 - 25 cm
34101-0012	Medium	25 - 27 cm
34101-0013	Large	27 - 30 cm
34101-0014	XLarge	30 - 33 cm



**Material:**


Beige/Black outer material:	Polyamid and Spandex
Foam:	Polyurethane
Edging:	Nylon & Spandex
Decoration:	Nylon & Spandex

**Washing Instruction:**

Close all Velcro prior to laundering.  
Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

**Contraindications:**

Edema, severe swelling.



# ELEMENTS


# BACK SUPPORT

## Instructions for ELEMENTS BACK SUPPORT Item no: 34102

**Intended use:** This product is intended for treatment of mild and temporary conditions such as pain, and a feeling of general weakness in the back. When wearing this support it gives some stabilization to the lower back due to compression. It also keeps the lower back warm.

### Sizes:

34102-0011	Small	82 - 95 cm
34102-0012	Medium	92 - 105 cm
34102-0013	Large	102 - 115 cm
34102-0014	XLarge	112 - 125 cm
34102-0015	XXLarge	122 - 135 cm



### Material:

Beige/Black outer material:	Polyamid & Spandex
Foam:	Polyurethane
Velcro:	Nylon
Stays in back:	POM
Stays on front:	Steel
Stay pockets:	Polyester
Eding:	Nylon & Spandex
Decroation:	Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering. Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

### Contraindications:

Back problems due to neurological origin.



# ELEMENTS


# WRIST SUPPORT

## Instructions for ELEMENTS WRIST SUPPORT Item No: 34103

**Intended use:** This product is intended for treatment of mild and temporary conditions such as pain, inflammation and joint stress in the wrist. The support limits the movement of the wrist, due to the aluminium stay on the palmar side. The support keeps the wrist warm.

### Sizes:

34103-1011	Small, left	11 - 14 cm
34103-1012	Medium, left	14 - 17 cm
34103-1013	Large, left	17 - 20 cm
34103-1014	XLarge, left	20 - 23 cm
34103-2011	Small, right	11 - 14 cm
34103-2012	Medium, right	14 - 17 cm
34103-2013	Large, right	17 - 20 cm
34103-2014	XLarge, right	20 - 23 cm



### Material:

Beige/Black outer material:	Polyamid & Spandex
Foam:	Polyurethane
Velcro:	Nylon
Support:	Aluminium
Support pocket:	Polyester
Edging:	Nylon & Spandex
Decoration:	Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering.  
Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

### Contraindications:

Severe edema, severe swelling.



# ELEMENTS


# EPI SUPPORT

## Instructions for ELEMENTS EPI SUPPORT Item No: 34104

**Intended use:** This product is intended for treatment of mild and temporary conditions of tennis elbow. It can also be combined with a wrist support for extended treatment. The support unloads the tension from the ligaments using light pressure over the muscle. The highest degree of function is achieved by using the two pads.

### Size:

34104-0000	Small/Medium	20 - 28 cm
34104-0013	Large/X-Large	28 - 36 cm



### Material:


Black outer material:	Polyamid & Spandex
Foam:	Polyurethane
In side material:	Polyester
Middle material:	Nylon
Velcro:	Nylon
D-ring:	POM
Edging:	Nylon & Spandex
Decoration:	Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering.  
Hand wash in luke warm water with a mild detergent. Use a washing bag. Dry flat. Do not tumble dry or dry clean.

### Contraindications:

Severe edema, severe swelling,



# ELEMENTS


# THUMB SUPPORT

## Instructions for ELEMENTS THUMB SUPPORT Item No: 34105

**Intended use:** This product is intended for treatment of mild and temporary conditions such as thumb pain, inflammation and joint stress. The support limits the movement of the thumb, and to some extent also the wrist. The support keeps the thumb and wrist warm.

### Sizes:

34105-1011	Small, left	< 17 cm
34105-1013	Large, left	> 17 cm
34105-2011	Small, right	< 17 cm
34105-2013	Large, right	> 17 cm



### Material:

Beige/Black outer material:	Polyamid & Spandex
Foam:	Polyurethane
Velcro:	Nylon
Plush:	Nylon & Polyurethane
Thumb strap:	Polyester
Edging:	Nylon & Spandex
Decoration:	Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering.  
Machine wash in 40°C (104°F). Use a washing bag. Dry flat. Do not tumble dry or dry clean.

### Contraindications:

Severe edema, severe swelling.



# ELEMENTS


# KNEE SUPPORT

## Instructions for ELEMENTS KNEE SUPPORT Item No: 34106

**Intended use:** This product is intended for treatment of mild and temporary conditions such as knee pain and inflammation. The support gives light compression and keeps the knee warm.

### Sizes:

34106-0010	XSmall	30 - 33 cm
34106-0011	Small	33 - 36 cm
34106-0012	Medium	36 - 39 cm
34106-0013	Large	39 - 42 cm
34106-0014	XLarge	42 - 47 cm



### Material:

Beige/Black outer material:	Polyamid & Spandex
Foam:	Polyurethane
Edging:	Nylon & Spandex
Decoration:	Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering.  
Machine wash in 40°C (104°F).  
Use a washing bag. Dry flat. Do not  
tumble dry or dry clean.

### Contraindications:

Edema, severe swelling, ligament injuries, leg  
ulcers.



# ELEMENTS

# BROMSA WRIST

## Instructions for ELEMENTS BROMSA WRIST Item No: 34107


**Intended use:** ELEMENTS BROMSA WRIST is intended to be used as prevention or treatment for lateral epicondylitis, overmobility, and unspecific pain in the wrist. This support is only intended for one person use.

Sizes: Small, Medium, Large, XLarge, left or right.

### Sizes:

34107-1011	Small, left	11 - 14 cm
34107-1012	Medium, left	14 - 17 cm
34107-1013	Large, left	17 - 20 cm
34107-1014	XLarge, left	20 - 23 cm

34107-2011	Small, right	11 - 14 cm
34107-2012	Medium, right	14 - 17 cm
34107-2013	Large, right	17 - 20 cm
34107-2014	XLarge, right	20 - 23 cm



### Material:

Beige/Black outer material: Polyamid & Spandex

Foam: Polyurethane

Velcro: Nylon

Support: Aluminium

Support pocket: Polyester

Edging: Nylon & Spandex

Decoration: Nylon & Spandex

### Washing Instruction:

Close all Velcro prior to laundering.

Remove plastic spoons.

Machine wash in 40°C (104°F).

Use a washing bag. Dry flat.

Do not tumble dry or dry clean.

### Contraindications:

Severe edema, severe swelling.



# ELEMENTS


# FOOT

## Instructions for ELEMENTS STRETCH UP Item No: 34108

**Intended use:** ELEMENTS STRETCH UP is a sock intended to be used when resting or at night. The sock is designed to keep tension on the tissue (plantar fasciitis) so it heals in a stretched position at night. It comes in three sizes, Small, Medium and Large. Fits both left and right.

**Sizes:**

34108-0011	Small	< 20 cm
34108-0012	Medium	20 - 24 cm
34108-0013	Large	> 40 cm



**Material specifications:**

This product is made from: 100% Spandex.

This product does not contain Latex, and is not flame retardant.

**Washing Instruction:**

Close all Velcro prior to laundering. Machine wash in 40°C (104°F).

Use a washing bag. Dry flat. Do not tumble dry or dry clean. Wash with similar colours.

**Contraindications:**

Circulation problems, open wounds on lower leg or foot, diabetes, skin problems or during pregnancy. Severe edema, Severe Swelling. Should not be used directly after any kind of surgery in leg or foot.



# ELEMENTS

# CERTIFICATE

## TESTEX AG, Swiss Textile Testing Institute

Gotthardstrasse 61  
CH-8027 Zürich

Institute of the International Association for Research  
and Testing in the Field of Textile Ecology



The company

### Allmay International Limited

RM 1405-06, 14/F.,  
Fu Lee Commercial Bldg.,  
14-20 Pilgrim Street, Jordan  
HK - Kowloon

is granted authorisation according to Oeko-Tex® Standard 100 to use the  
Oeko-Tex® mark, based on test report HKA0 085398.1



**Tested for harmful substances**  
according to Oeko-Tex® Standard 100  
HKA0 060904 TESTEX

for the following articles:

**Grey polyurethane foam laminated with knitted  
polyamide/spandex fabric, dyed (for article No. 34100 Elements  
ANKLE/ 34101 Elements ELBOW/ 34102 Elements BACK/ 34103  
Elements WRIST/ 34104 Elements TENNIS/ 34105 Elements  
THUMB/ 34106 Elements KNEE) (based on materials partly pre-  
certified according to Oeko-Tex® Standard 100)**

The results of the inspection made according to Oeko-Tex® Standard 100, product class II for  
products with direct contact to skin, have shown that the above mentioned goods meet the  
human-ecological requirements of the standard presently established.

The certified articles fulfil the requirements of Annex XVII of REACH (incl. the use of azo-  
dyes, nickel, etc.) as well as the American requirement regarding total content of lead in  
children's articles (CPSIA; with the exception of accessories made from glass).

The holder of the certificate, who has issued a conformity declaration according to EN ISO  
17050-1, is obliged to use the Oeko-Tex® mark only in conjunction with products that are  
conform with the sample initially tested.

**This authorisation is valid until 15.11.2013**

Zurich, 01.10.2012



*Serge Rolle*

Serge Rolle  
CEO

*Jean Pierre Haug*

Dr. Jean Pierre Haug  
COO





# *ELEMENTS*

**SVERIGE**

Tel 042-25 27 01  
Fax 042-25 27 25  
mail.sweden@camp.se

**DANMARK**

43 96 66 99  
43 43 22 66  
mail.denmark@camp.dk

**FINLAND**

09-350 76 30  
09-350 76 338  
mail.finland@camp.fi

**NORGE**

23 23 31 20  
23 23 31 21  
mail.norway@camp.no