

VISSION™ Neo Thumb Long

Soft and compliant thumb orthosis in latex-free Neoprene that provides both warmth and support to CMC and MCP joints. At rest, it keeps the thumb in a neutral position. Available in two lengths. The long model provides support for the wrist as well. The material can be cut without fraying.

Recommended Range Of Application

CMC 1 osteoarthritis, rheumatoid arthritis, after minor trauma, in case of hypermobility in the thumb and wrist, after arthroplasty surgery and in case of overload.

Measurements

Wrist circumference



Item No.	Model	Size	L/R	Color
460261000	Long	Onesize	Left	Black
460262000	Long	Onesize	Right	Black