

VISSION™ Thumb Orthosis Short with Stay

This thumb orthosis have a plastic stay on the thumb side to support the thumb and the surrounding structures. The material under the stay gives good padding to avoid pressure on the joint or muscles. The material is durable and stable, with a comfortable inside. The thumb strap is open for adjustability.

Recommended Range Of Application

This product is intended to be used for instability in the thumb joints (CMC- or MCP 1 joints) and to position the thumb in a functional positions. Indications such as inflammation, joint stress and for rheumatic degenerations.



Item No.	Model	Size	Measurement	L/R	Color
460281009	Short with stay	XX-Small	11 - 12 cm	Left	Black
460281010	Short with stay	X-Small	12 - 13 cm	Left	Black
460281011	Short with stay	Small	13 - 14 cm	Left	Black
460281012	Short with stay	Medium	14 - 16 cm	Left	Black
460281013	Short with stay	Large	16 - 19 cm	Left	Black
460281014	Short with stay	X-Large	19 - 22 cm	Left	Black
460282009	Short with stay	XX-Small	11 - 12 cm	Right	Black
460282010	Short with stay	X-Small	12 - 13 cm	Right	Black
460282011	Short with stay	Small	13 - 14 cm	Right	Black
460282012	Short with stay	Medium	14 - 16 cm	Right	Black
460282013	Short with stay	Large	16 - 19 cm	Right	Black
460282014	Short with stay	X-Large	19 - 22 cm	Right	Black