

# BASIC Wrist with Thumb

Stable Wrist support that also supports the thumb without too much material between the thumb and the fingers. The thumb section is partly free from the other material to allow for a good support around the thumb structures. The volar aluminum splint can be removed. A narrow plastic stay run along the side of the thumb for extra support.

## Recommended Range Of Application

Post-traume, Rheumatic problems, instability and pain releif.

## Measurements

Wrist circumference



Item No.	Size	Measurement	L/R	Color
274421010	X-Small	< 12 cm	Left	Blue
274421011	Small	11-14 cm	Left	Blue
274421012	Medium	14-17 cm	Left	Blue
274421013	Large	17-20 cm	Left	Blue
274421014	X-Large	20-23 cm	Left	Blue
274422010	X-Small	< 12 cm	Right	Blue
274422011	Small	11-14 cm	Right	Blue
274422012	Medium	14-17 cm	Right	Blue
274422013	Large	17-20 cm	Right	Blue
274422014	X-Large	20-23 cm	Right	Blue