

VISSION™ Shoulder Orthosis Immo

VISSION™ Shoulder Orthosis Immo is designed to immobilize and support the upper limb, with the elbow flexed at 90° or 130° as needed. It is particularly suited for shoulder and pectoral girdle injuries. The product is reversible for use on either shoulder and intended for single-patient use.

Recommended Range Of Application

Shoulder, arm, elbow and pectoral girdle trauma (fractures, dislocations and sprains). Shoulder and pectoral girdle surgery. Neurology: paralysis of brachial-plexus.

Contraindications

Non known.

Measurements

Measure 1 = Chest circumference
Measure 2 = Length upper arm



Item No.	Model	Size	Measurement 1	Measurement 2
460220010	Shoulder Orthosis IMMO	X-Small	68 - 78 cm	< 30 cm
460220011	Shoulder Orthosis IMMO	Small	80 - 90 cm	< 38 cm
460220012	Shoulder Orthosis IMMO	Medium	90 - 100 cm	< 38 cm
460220013	Shoulder Orthosis IMMO	Large	100 - 113 cm	< 38 cm
460220014	Shoulder Orthosis IMMO	X-Large	113 - 120 cm	> 38 cm
460220015	Shoulder Orthosis IMMO	XX-Large	120 - 140 cm	> 38 cm