

# Vission™ Shoulder Orthosis IMMO

Vission shoulder immobilizer is intended to immobilize and contain the upper limb, elbow, forearm flexed at 90 or 130 degrees as required. Especially suited to shoulder and pectoral girdle lesions. The product is reversible and can be used either on the left or on the right shoulder as required. The product is intended for one person use.

### Recommended Range Of Application

Shoulder, arm, elbow and pectoral girdle trauma (fractures, dislocations and sprains). Shoulder and pectoral girdle surgery. Neurology: paralysis of brachial-plexus.

### Contraindications

Non known.

### Measurements

Measure 1 = Chest circumference

Measure 2 = Length upper arm



Item No.	Model	Size	Measurement 1	Measurement 2
460220010	Shoulder Orthosis IMMO	X-Small	68 - 78 cm	< 30 cm
460220011	Shoulder Orthosis IMMO	Small	80 - 90 cm	< 38 cm
460220012	Shoulder Orthosis IMMO	Medium	90 - 100 cm	< 38 cm
460220013	Shoulder Orthosis IMMO	Large	100 - 113 cm	< 38 cm
460220014	Shoulder Orthosis IMMO	X-Large	113 - 120 cm	> 38 cm
460220015	Shoulder Orthosis IMMO	XX-Large	120 - 140 cm	> 38 cm