

# COMBO™

Dynamic low profile, lightweight, functional orthotic solution for the management of genu recurvatum or chronic knee instability, accompanied by footdrop. Low profile and lightweight – a KAFO your patients will wear! Uprights easily shaped and adjusted to optimum height for patient. Modular components allow you to accommodate variances in thigh and calf circumference. Easy to don and doff. Easily removed when wearer is relaxing (watching TV, etc.)



## Recommended Range Of Application

Knee instability, knee hyperextension in combination with a foot drop from stroke, post-polio, MS, etc.

## Contraindications

Use of COMBO with any product other than ToeOFF® or BlueROCKER® is prohibited.

## Measurements

Kneemeasure 1: M-L mid-patella

Kneemeasure 2: Circumference 20 cm above the mid-patella

Kneemeasure 3: M-L 20 cm above the mid-patella

| Item No.  | Model                                     | Size     | L/R       | Knee Measurement 1 | Knee Measurement 2 | Knee Measurement 3 |
|-----------|---|----------|-----------|--------------------|--------------------|--------------------|
| 287800000 | Joints with uprights                      | One size | Bilateral | -                  | -                  | -                  |
| 287811011 | Thigh cuff                                | Small    | Left      | -                  | 30-40 cm           | 11-15 cm           |
| 287811012 | Thigh cuff                                | Medium   | Left      | -                  | 38-50 cm           | 13-17 cm           |
| 287811013 | Thigh cuff                                | Large    | Left      | -                  | 48-60 cm           | 15-19 cm           |
| 287811014 | Thigh cuff                                | X-Large  | Left      | -                  | 58-70 cm           | 17-21 cm           |
| 287812011 | Thigh cuff                                | Small    | Right     | -                  | 30-40 cm           | 11-15 cm           |
| 287812012 | Thigh cuff                                | Medium   | Right     | -                  | 38-50 cm           | 13-17 cm           |
| 287812013 | Thigh cuff                                | Large    | Right     | -                  | 48-60 cm           | 15-19 cm           |
| 287812014 | Thigh cuff                                | X-Large  | Right     | -                  | 58-70 cm           | 17-21 cm           |
| 287820011 | Strap set and buckles, poplitea interface | Small    | Bilateral | 7-10,5 cm          | -                  | -                  |
| 287820012 | Strap set and buckles, poplitea interface | Medium   | Bilateral | 9,5-13 cm          | -                  | -                  |
| 287820013 | Strap set and buckles, poplitea interface | Large    | Bilateral | 12,5-16 cm         | -                  | -                  |
| 287960000 | Calf Strap                                | One Size | -         | -                  | -                  | -                  |
| 287830011 | Thigh cuff softKIT, 2 pc.                 | Small    | Bilateral | -                  | -                  | -                  |
| 287830012 | Thigh cuff softKIT, 2 pc.                 | Medium   | Bilateral | -                  | -                  | -                  |
| 287830013 | Thigh cuff softKIT, 2 pc.                 | Large    | Bilateral | -                  | -                  | -                  |
| 287830014 | Thigh cuff softKIT, 2 pc.                 | X-Large  | Bilateral | -                  | -                  | -                  |
| 287840011 | Tibia connector                           | Small    | Bilateral | -                  | -                  | -                  |
| 287840012 | Tibia connector                           | Medium   | Bilateral | -                  | -                  | -                  |
| 287840013 | Tibia connector                           | L-XL     | Bilateral | -                  | -                  | -                  |