

# Vission™ Thumb Orthosis Long with Stay

Soft and elastic wrist orthosis in a long model. The orthosis has a removable plastic stay that gives support to the thumb and surrounding area as well as light support for the wrist. This thumb orthosis is intended to restrict movements in CMC and MCP 1 joints and to position the thumb in a functional position. It is made of a durable and resilient material with a soft and comfortable lining. Open adjustable thumb grip.

### Recommended Range Of Application

CMC and MCP 1 osteoarthritis or instability, Trigger Finger, Morbus De Quervain.



Item No.	Model	Size	Measurement	L/R	Color
460231009	Long with stay	XX-Small	11 - 12 cm	Left	Black
460231010	Long with stay	X-Small	12 - 13 cm	Left	Black
460231011	Long with stay	Small	13 - 14 cm	Left	Black
460231012	Long with stay	Medium	14 - 16 cm	Left	Black
460231013	Long with stay	Large	16 - 19 cm	Left	Black
460231014	Long with stay	X-Large	19 - 22 cm	Left	Black
460232009	Long with stay	XX-Small	11 - 12 cm	Right	Black
460232010	Long with stay	X-Small	12 - 13 cm	Right	Black
460232011	Long with stay	Small	13 - 14 cm	Right	Black
460232012	Long with stay	Medium	14 - 16 cm	Right	Black
460232013	Long with stay	Large	16 - 19 cm	Right	Black
460232014	Long with stay	X-Large	19 - 22 cm	Right	Black