

## ToeOFF® Flow 2 ½

Using our experience and development in manufacturing technology, together with your feedback, we have developed a new Allard AFO - ToeOFF\$ Flow 2  $\frac{1}{2}$ .

ToeOFF® Flow 2  $\frac{1}{2}$  has a unique formula with new features and benefits which include:

½ heel height of 7mm in all sizes, more ROM in A-P plane and smoother surface underneath the footplate. The low profile tibia plate is more dynamic in the A-P plane and it has a better fit in today's lower heeled shoes.

ToeOFF Flow  $2\frac{1}{2}$  has all Microfix pre-applied which makes it quicker to attach a SoftKIT  $^{\text{TM}}$  2.0, ComfortKIT  $^{\text{TM}}$  2.0 or CoverKIT  $^{\text{TM}}$  2.0. Closure straps are the same as for ToeOFF $^{\text{R}}$ 2.0

## **Recommended Range Of Application**

Foot drop combined with moderate to no spasticity. Limb proprioception deficit and mild proximal deficit. For further information on uses and contraindications see AFO Professional Instruction.

## **Contraindications**

Severe edema, leg ulcers, moderate to severe spasticity, medial-lateral instability.







| Item No.  | Model                         | Size    | L/R   | Footplate Length | Height  |
|-----------|-------------------------------|---------|-------|------------------|---------|
| 289601010 | ToeOFF® Flow2½                | X-Small | Left  | 21 cm            | 36 cm   |
| 289601011 | ToeOFF® Flow2½                | Small   | Left  | 23 cm            | 38 cm   |
| 289601012 | ToeOFF® Flow2½                | Medium  | Left  | 24,5 cm          | 40,5 cm |
| 289601013 | ToeOFF® Flow2½                | Large   | Left  | 27 cm            | 43 cm   |
| 289601014 | ToeOFF® Flow2½                | X-Large | Left  | 28,5 cm          | 45 cm   |
| 289602010 | ToeOFF® Flow2½                | X-Small | Right | 21 cm            | 36 cm   |
| 289602011 | ToeOFF® Flow2½                | Small   | Right | 23 cm            | 38 cm   |
| 289602012 | ToeOFF® Flow2½                | Medium  | Right | 24,5 cm          | 40,5 cm |
| 289602013 | ToeOFF® Flow21/2              | Large   | Right | 27 cm            | 43 cm   |
| 289602014 | ToeOFF® Flow2½                | X-Large | Right | 28,5 cm          | 45 cm   |
| 289621010 | ToeOFF® Flow 2 ½ with D-Rings | X-Small | Left  | 21 cm            | 36 cm   |
| 289621011 | ToeOFF® Flow 2 ½ with D-Rings | Small   | Left  | 23 cm            | 38 cm   |
| 289621012 | ToeOFF® Flow 2 ½ with D-Rings | Medium  | Left  | 24,5 cm          | 40,5 cm |
| 289621013 | ToeOFF® Flow 2 ½ with D-Rings | Large   | Left  | 27 cm            | 43 cm   |
| 289621014 | ToeOFF® Flow 2 ½ with D-Rings | X-Large | Left  | 28,5 cm          | 45 cm   |
| 289622010 | ToeOFF® Flow 2 ½ with D-Rings | X-Small | Right | 21 cm            | 36 cm   |
| 289622011 | ToeOFF® Flow 2 ½ with D-Rings | Small   | Right | 23 cm            | 38 cm   |
| 289622012 | ToeOFF® Flow 2 ½ with D-Rings | Medium  | Right | 24,5 cm          | 40,5 cm |
| 289622013 | ToeOFF® Flow 2 ½ with D-Rings | Large   | Right | 27 cm            | 43 cm   |
| 289622014 | ToeOFF® Flow 2 ½ with D-Rings | X-Large | Right | 28,5 cm          | 45 cm   |