

allard_{NT}

Humerux

offers maximum control thanks to the innovative traction system.

traction
system®



Humerux

offers maximum control thanks to the innovative traction system.

Indicated for the treatment of diaphyseal humeral fractures, postoperative treatment of fractures and control of the instability of the fractures focus during the consolidation process.




▲ **Y-traction System uses a traction device to increase and adjust the tension**, helping to maximise control over the alignment, fixation and immobilisation of the humerus



Ytraction system®

▲ **Padded and breathable technical fabric with comfort and smoothness towards the skin!**

100% LATEX FREE ✓



▲ The plates and straps stabilize the humerus by applying uniform and personalized compression to the surrounding soft tissue.

The bi-valve structure in line with the proven Sarmiento method gives equal compression around the soft tissue. These can be individually shaped*.



▲ Easy-fit Chest Straps to prevent the orthosis from sliding down. Padded axilla support for comfortable underarm fitting, avoiding any type of friction or discomfort.

▲ Fixation system using 3 humeral straps, threaded through 3 butterfly buckles with micro-hook closure and an anchor stop to ensure an easy, comfortable and secure adjustment.

*Remove the thermoplastic plates from the textile. They can then be heated with a heatgun at 150°C with about 15 cm distance.

Item No.	Size	Length	Circumference
62105 0011	Small	37 cm	18 - 26 cm
62105 0012	Medium	38 cm	26 - 34 cm
62105 0013	Large	41 cm	34 - 42 cm
62105 0014	X-Large	42 cm	42 - 50 cm

