

# SELECTION® Open Wrist Children

SELECTION® Open Wrist Children is a brace that is simple to put on and take off. It can be opened totally flat to make it easier to put on. The integrated "wrap around" straps gives increased stability around the wrist. It has two thumb straps. One around the thumb and another on the dorsal side of the hand. This keeps the thumb in a correct position so that it does not drop towards the palm of the hand.

The orthosis is intended to immobilize the CMC joint and wrist in order to reduce pain, to guide the wrist into a neutral position and to support a wrist after trauma/post operatively if needed.

- The two-piece thumb strap keeps the thumb in a neutral and functional position.
- Integrated stabilizing strap around the wrist.
- The extra material supports the thumbs muscles and the extra strap helps to move the adducted thumb out from the palm.



## Recommended Range Of Application

Rheumatoid Arthritis, CP.

Item no. Left	Item no. Right	Size	Wrist circumference	Palm length
353061011	363062011	Small	11,5- 14 cm	7-8 cm
363061012	363062012	Medium	12 - 15,5 cm	8-9 cm
363061013	363062013	Large	14 - 16,5 cm	8-9 cm