

# SELECTION® Wrist Rigid

Selection Wrist Rigid has the same features as the soft version, but also a thumb section added for increased thumb support and stability. The thumb strap can be individually adjusted. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. This product is intended to be used when the first CMC-joint needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear.

## Recommended Range Of Application

Inflammations, work related strain or when rheumatic and/or degenerative changes are present.

## Contraindications

Non known.

## Measurements

Wrist circumference

SELECTION®



| Item No.  | Model | Size     | Measurement | L/R   | Color |
|-----------|-------|----------|-------------|-------|-------|
| 352061010 | Rigid | X-Small  | < 11 cm     | Left  | Black |
| 352061011 | Rigid | Small    | 11 - 14 cm  | Left  | Black |
| 352061012 | Rigid | Medium   | 14 - 17 cm  | Left  | Black |
| 352061013 | Rigid | Large    | 17 - 20 cm  | Left  | Black |
| 352061014 | Rigid | X-Large  | 20 - 23 cm  | Left  | Black |
| 352061015 | Rigid | XX-Large | 23 - 26 cm  | Left  | Black |
| 352062010 | Rigid | X-Small  | < 11 cm     | Right | Black |
| 352062011 | Rigid | Small    | 11 - 14 cm  | Right | Black |
| 352062012 | Rigid | Medium   | 14 - 17 cm  | Right | Black |
| 352062013 | Rigid | Large    | 17 - 20 cm  | Right | Black |
| 352062014 | Rigid | X-Large  | 20 - 23 cm  | Right | Black |
| 352062015 | Rigid | XX-Large | 23 - 26 cm  | Right | Black |