

ToeOFF® 2.0

ToeOFF® 2.0 is the updated version of the Original ToeOFF®. This new model has several new features that have been asked for from customers over time. The shape is more organic, with shorter side extensions and more rounded shape. Also the surface is smoother and more scratch resistant. The MikroFix is pre-attached on the inside of the brace to accomodate a faster fitting procedure. It is also available with two different strap systems - With D-ring or without D-ring. ToeOFF® 2.0 has a heel height between 12 - 15 mm depending on size.

Recommended Range Of Application

ToeOFF® 2.0 is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy. ToeOFF® 2.0 is designed to support gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations. ToeOFF® 2.0 can be used for toe amputations.

Contraindications

ToeOFF® 2.0 should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

ToeOFF® 2.0



Item No.	Model	Size	L/R	Shoe Size	Footplate Length	Height
289201010	ToeOFF 2.0 Simple Strap	X-Small	Left	32-35	21 cm	36 cm
289201011	ToeOFF 2.0 Simple Strap	Small	Left	35-38	23 cm	38 cm
289201012	ToeOFF 2.0 Simple Strap	Medium	Left	38-42	24,5 cm	40,5 cm
289201013	ToeOFF 2.0 Simple Strap	Large	Left	42-45	27 cm	43 cm
289201014	ToeOFF 2.0 Simple Strap	X-Large	Left	>45	28,5 cm	45 cm
289202010	ToeOFF 2.0 Simple Strap	X-Small	Right	32-35	21 cm	36 cm
289202011	ToeOFF 2.0 Simple Strap	Small	Right	35-38	23 cm	38 cm
289202012	ToeOFF 2.0 Simple Strap	Medium	Right	38-42	24,5 cm	40,5 cm
289202013	ToeOFF 2.0 Simple Strap	Large	Right	42-45	27 cm	43 cm
289202014	ToeOFF 2.0 Simple Strap	X-Large	Right	>45	28,5 cm	45 cm
289221010	ToeOFF 2.0 D-Ring Strap	X-Small	Left	32-35	21 cm	36 cm
289221011	ToeOFF 2.0 D-Ring Strap	Small	Left	35-38	23 cm	38 cm
289221012	ToeOFF 2.0 D-Ring Strap	Medium	Left	38-42	24,5 cm	40,5 cm
289221013	ToeOFF 2.0 D-Ring Strap	Large	Left	42-45	27 cm	43 cm
289221014	ToeOFF 2.0 D-Ring Strap	X-Large	Left	>45	28,5 cm	45 cm
289222010	ToeOFF 2.0 D-Ring Strap	X-Small	Right	32-35	21 cm	36 cm
289222011	ToeOFF 2.0 D-Ring Strap	Small	Right	35-38	23 cm	38 cm
289222012	ToeOFF 2.0 D-Ring Strap	Medium	Right	38-42	24,5 cm	40,5 cm
289222013	ToeOFF 2.0 D-Ring Strap	Large	Right	42-45	27 cm	43 cm
289222014	ToeOFF 2.0 D-Ring Strap	X-Large	Right	>45	28,5 cm	45 cm