VISSION™ Finger Splint Multi



For immobilizing the fingers following trauma or operation and preventing contractures of the finger joints. Can also be used in conservative finger therapy and the treatment of Dupuytren syndrome. Conditions requiring temporary finger immobilization. The splint can also be used in connection with the treatment of fractures on metacarpus III-V and its associated phalanges.

Recommended Range Of Application
Immobilizing, Dupuytren syndrome and fractures.

Contraindications
Non known.

Measurements
Wrist circumference.

Item No.	Model	Size	Measurement	L/R	Color
460211011	Vission Multifinger	Small	< 16 cm	Left	Black
460211012	Vission Multifinger	Medium	16 - 19 cm	Left	Black
460211013	Vission Multifinger	Large	> 19 cm	Left	Black
460212011	Vission Multifinger	Small	< 16 cm	Right	Black
460212012	Vission Multifinger	Medium	16 - 19 cm	Right	Black
460212013	Vission Multifinger	Large	> 19 cm	Right	Black