# VISSION™ Shoulder Orthosis Immo



VISSION ™ Shoulder Orthosis Immo is designed to immobilize and support the upper limb, with the elbow flexed at 90° or 130° as needed. It is particularly suited for shoulder and pectoral girdle injuries. The product is reversible for use on either shoulder and intended for single-patient use.

## Recommended Range Of Application

Shoulder, arm, elbow and pectoral girdle trauma (fractures, dislocations and sprains). Shoulder and pectoral girdle surgery. Neurology: paralysis of brachial-plexus.

## Contraindications

Non known.

## Measurements

Measure 1 = Chest circumference
Measure 2 = Length upper arm

| Item No. | Size | Measurement 1 | Measurement 2 |
|---|---|---|---|
| 460220010 | X-Small | 68 - 78 cm | < 30 cm |
| 460220011 | Small | 80 - 90 cm | < 38 cm |
| 460220012 | Medium | 90 - 100 cm | < 38 cm |
| 460220013 | Large | 100 - 113 cm | < 38 cm |
| 460220014 | X-Large | 113 - 120 cm | > 38 cm |
| 460220015 | XX-Large | 120 - 140 cm | > 38 cm |