SELECTION® Radialis Thumb Part



Supporting textile thumb section to be combined with SELECTION® Wrist orthoses 35204, 35205, 35209, 35222, 35223, or 35305. Read the instructions for the main product before use. In case of radial nerve palsy (wrist drop caused by damage to the radial nerve), the thumb section can be combined with the 35217 SELECTION® Radial Finger Part to enable a functional grip while the nerve recovers.

Available in smaller sizes for children
The Thumb Part is also available in XS and XXS, making it compatible with our children's wrist orthosis. This makes it particularly well-suited for children with conditions such as cerebral palsy, who often require more tailored solutions.

Recommended Range Of Application

For relief and guidance of the thumb, Radial Nerve Palsy.

Measurements

Length joint centre MCP1 – joint centre IP1

Item No.	Model	Size	Measurement	L/R	Color
352161009	Radialis Thumb Part	XX-Small	1,5-2 cm	Left	Black
352161010	Radialis Thumb Part	X-Small	2-3 cm	Left	Black
352161011	Radialis Thumb Part	Small	3 - 4 cm	Left	Black
352161013	Radialis Thumb Part	Large	5 - 6 cm	Left	Black
352162009	Radialis Thumb Part	XX-Small	1,5-2 cm	Right	Svart
352162010	Radialis Thumb Part	X-Small	2-3 cm	Right	Black
352162011	Radialis Thumb Part	Small	3 - 4 cm	Right	Black
352162013	Radialis Thumb Part	Large	5 - 6 cm	Right	Black